## ENGLISH

### INSTRUCTIONS FOR USE:

This package insert is not a reference to surgical techniques. It is designed to assist in using this product.

### PRODUCT DESCRIPTION:

LUXCRYL MONOFAST is a synthetic absorbable sterile surgical suture composed of Poliglecaprone 25, a copolymer of glycolide and  $\epsilon$ -Caprolactone (poly(glycolide-co-caprolactone)). The monofilament sutures are undyed, in the natural beige colour. Sutures are also available dyed in violet.

Progressive loss of tensile strength, mass and eventual absorption of the sutures occur by means of hydrolysis. The loss of effective tensile strength is approximately 25% 1 week post-implantation. The absorption is complete after 90 to 110 days. LUXCRYL MONOFAST is available in a range of gauge sizes and lengths, non-needed or attached to atraumatic stainless steel needles of varying types. The sutures are individually packed in sterile foil packs.

LUXCRYL MONOFAST complies the requirements of the European Pharmacopoeia Monograph and the United States Pharmacopoeia.

**INDICATIONS:** LUXCRYL MONOFAST is indicated for use in general soft tissue approximation.

The selection of the suture should be based on the patient's condition, the surgical experience, the surgical technique to be employed and the wound to be treated.

### CONTRA-INDICATIONS:

LUXCRYL MONOFAST should not be used:

- in tissues which may undergo expansion, stretching or distension or which may require a long term mechanical support.
- in cardiovascular and neurological tissues,
- in microsurgery and in ophthalmic surgery,
- in the closure of fascial tissue,
- in patients that are allergic to suture material.

The use of this suture can be inappropriate in elders, patients that are not properly fed, weak patients and patients under certain conditions that could delay the healing of the wound.

### ADVERSE REACTIONS:

Adverse reactions associated with this device include: wound dehiscence, dehiscence of the plate, failure to furnish a sutured or distended wound, healing, transient irritation at the wound site, transient inflammatory foreign body response, septic infection, erythema and the wound side and calculi formation.

Like all foreign bodies LUXCRYL MONOFAST may potentiate an existing infection.

### HANDLING INSTRUCTIONS / WARNTINGS / PRECAUTIONS:

Users should be familiar with surgical procedures and techniques and the techniques of suture implying absorbable sutures before employing LUXCRYL MONOFAST suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.

Particular factors of each patient should be considered regarding healing process in vivo.

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Les sutures cutanées, vaginales et conjonctivales qui doivent être snipped off ou removed as indicated. Place subcuticular sutures as deeply as possible to minimize the erythema and induration. Avoid prolonged contact of any suture with salt solutions (urinary or biliary), as it may result in calculus formation. Contaminated or infected wounds should be managed with acceptable surgical practice. A supplemental nonabsorbable suture or an immobilisation of joints by an external support should be considered by the surgeon in the closure of the sites which may require additional support (e.g. expansion, distension, etc.). In tissues with poor blood supply extrusion and delayed absorption may occur.

Open the outer bag of the foil pack by pulling its two loose ends away from each other – in such a way that the sterile bag containing the suture remains intact until the opening of the sterile surface. The sterile bag should then only be touched with sterile gloves and opened by tearing the upper part of the bag.

The product is intended for single use!

When suturing, remember:

• Always double check that you have tied a tight knot. Additional throws should be considered in some surgical circumstances.

• Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury, which may for example result in transmission of blood-borne pathogens.

• When suturing be careful to avoid damage (suture damage, needle damage, etc.) with surgical instruments (e.g. needle holder, forceps, etc.) or needle.

• Hold the needle in the area between the third and the half of the distance between the extreme of the union with the thread and the tip. If you hold it near to the tip, you could damage the functional integrity of the needle or break it. If you hold it by the union zone with the thread, the needle body could be broken. When the needle function is altered, it could generate a decrease in its folding resistance and break.

• Discard any unused suture material remaining. Suture needles should be discarded into separated "sharps" container. Dispose of contaminated devices and packaging materials utilizing standard hospital procedures and universal precautions for bio-hazardous waste.

• Do not use rusty needles! Use of rusty needles may cause infection and other harms.

Do not use if the package is opened or damaged.

This device was designed, tested and manufactured for single patient use only. Reuse, use of the device with opened or damaged packaging, reprocessing and/or resterilization of this device may lead to its failure and subsequent injury, illness or death of the patient or render the risk of contamination and patient infection, illness or death of the patient.

### STERILISATION:

LUXCRYL MONOFAST is sterilised by ethylene oxide.

**STORAGE:** The product should be protected from direct sunlight and heat, and stored in its original packaging in a clean, dry room at a temperature from 5°C to 30°C.

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