ENGLISH
INSTRUCTIONS FOR USE:
This package insert is not a reference to surgical techniques. It is designed to assist in using this product.

PRODUCT DESCRIPTION:

LUXCRYL PGA RAPID is a synthetic absorbable sterile surgical suture composed of polyglycolic acid. The sutures are coated with a hydrophilic polypropylene and calcium polyglycolate. Ces fils braided sutures are moist and dyed in violet. Progressive loss of tensile strength, mass and eventual absorption of the sutures occur by means of hydrolysis. The loss of effective tensile strength is approximately 50% after 7 days post-implantation. The absorption is complete after approximately 42 days.

LUXCRYL PGA RAPID is available in a range of gauge sizes and lengths, non-needled or attached to atraumatic stainless steel or titanium carrying types. The sutures are individually packed in sterile foil packs.

LUXCRYL PGA RAPID complies the requirements of the European Pharmacopoeia Monograph and the United States Pharmacopoeia.

INDICATIONS:
LUXCRYL PGA RAPID est conforme à la monographie de la LUXCRYL PGA RAPID is indicated only for use in superficial lacerations and wounds on the arm and hand support is required and where the absorption of the suture would be beneficial. Due to its absorption profile LUXCRYL PGA RAPID is useful for skin closure.

The selection of the suture should be based on the patient's condition, the surgical experience, the surgical technique to be employed and the wound to be treated.

CONTRA-INDICATIONS:

LUXCRYL PGA RAPID should not be used:

- where extended approximation of tissues under stress is required,
- when wound support beyond 7 days is required,
- for ligation,

- in tissues which may undergo expansion, stretching or distension such as the following:
 - in cardiovascular and neurological tissues,
 - in microsurgery and ophthalmic surgery,
 - in patients that are allergic to suture material.
- The use of this suture can be improper in elders, patients that are not properly fed, weak patients, and patients under certain conditions that could delay the healing of the wound.

ADVERSE REACTIONS associated with this device include: wound dehiscence, failure to provide adequate wound support in sites where expansion, stretching or distension occur or in patients with delayed wound healing, transient local irritation at the wound site, transient inflammatory foreign body response, induration, erythema at the wound site and calculi formation. Use of the suture bodies LUXCRYL PGA RAPID may potentiate an existing infection.

HANDLING INSTRUCTIONS / WARNINGS / PRECAUTIONS: Users should be familiar with surgical procedures and techniques and tie techniques involving absorbable sutures before employing LUXCRYL PGA RAPID suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.

Patient and patient's caregiver should be considered regarding- healing process in vivo.

Skin, vaginal and conjunctival sutures which must remain in place longer than 7 days may cause localised irritation and should be snipped off or removed as indicated. Place subcuticular sutures as deep as possible to minimize the erythema and induration. Avoid prolonged contact of any suture with all solutions (urine, sweat, etc.) as this may result in a localised infection. Contaminated or infected wounds should be managed with acceptable surgical practice. A supplemental nonabsorbable suture or an immobilisation of joints by an external support should be considered by the surgeon in the closure of the sites which may require additional support (e.g. expansion, distension, etc.).

In tissues with poor blood supply suture extrusion and delayed wound healing may occur. Patient should be considered regarding- healing process in vivo.

Always double check that you have tied a tight knot. Additional threads should be considered in some surgical circumstances.

Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury, which may for example result in transmission of blood-borne pathogens.

When suturing be careful to avoid damage (suture damage, needle damage, etc) with surgical instruments (e.g. needle holder, forceps, etc.) or needles.

Hold the needle in the area between the third and the half of the distance between the extreme of the union with the thread and the tip. If you hold it near to the tip, you could damage the functional integrity of the needle or break it. If you hold it by the union zone with the thread, the needle body could be broken.

When the needle function is altered, it could generate a decrease in its folding resistance and breaking strength.

Discard any unused suture material remaining. Suture needles should be discarded into separated "sharps" container. Dispose of contaminated devices and packaging materials utilizing standard hospital procedures and universal precautions for bio hazardous waste.

Do not use "sticky needles! Use of unused needles may cause infection and other harms.

Do not restitute!

This device was designed, tested and manufactured for single patient use only. Reuse, use of the device with opened or damaged packaging, reprocessing and/or restenilization of this device may lead to its failure and subsequent injury, illness or death of the patient and create the risk of contamination and patient infection, illness or death of the patient.

STERILISATION:
LUXCRYL PGA RAPID is sterilised by ethylene oxide.

STORAGE:
The product should be protected from direct sunlight and heat, and stored in its original packaging in a clean, dry room at a constant temperature from 5°C to 30°C.

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